


REGIONAL ANAESTHESIA



**26th Year
Manufacturing
Since 1993**





M.D.L. Srl
Via Tavani 1/a 23014 Delebio (So) Italy
Tel +39 0342 682130 Fax +39 0342 691316
email : info@mdlsrl.com web : www.mdlsrl.com




INDEX

STIMYPLUS	06
STIMYULTRA	07
STIMYECHO	08
STIMY	09
P-RF	10
Tables	11

StimyPlus


Insulated nerve block needle for
combined single shot E.N.S.
Technique and Ultrasound Technique



StimyPlus is a hybrid needle with the ability to be used with both ENS and ultrasound techniques.

Electrical wire with universal 2 mm connector

Detachable luer extension tube



pp.12

REF: AETMxxyyy

16°



REF: AEQMxxyyy

30°



StimyUltra

Insulated nerve block needle,
single shot E.N.S. technique



pp.11

REF: AETRxxxxyy 16°


REF: AEQRxxxxyy 30°

REF: AEYRxxxxyy 20°



StimyEcho


*Insulated nerve block needles,
single shot Ultrasound Technique*



Nanoclear coated
cannula with centimeter
depth marking

StimyEcho is for single shot
nerve blocks using only
ultrasound technique.

Detachable luer
extension tube



pp.12

REF: AETMNxxyyy

16°

REF: AEQMNxxyyy

30°

pp.13



Stimy

*Insulated nerve block needle,
single shot E.N.S. technique*

2.4" LED backlit
color **LCD screen**

Shortcut buttons
allow direct switching
of key parameter
values during
procedure

**Rotating
knob** allows
precise setting
of operating
parameters

Ergonomic design allows
for easy handling

**Acoustic signals warn
when parameters are out of
acceptance limit**



P-RF

Radiofrequency introducer cannulas
for Pain Therapy

REF: RF

pp.13/14

Indicated for percutaneous introduction of probes in the procedure of radiofrequency pulsed or lesive in articular and intraarticular facets. The needle without coating is indicated for percutaneous drug's injection in articular and intraarticular facets.



Indicated for percutaneous introduction of probes in the procedure of intraforaminal radiofrequency pulsed or lesive. The needle without coating is indicated for percutaneous drug's injection intraforaminal.



Catheter introduction with special spoon point in order to prevent any damage during introduction and extraction



ORDER INFORMATION

Stimy Ultra 16 - pp.07

CODE	Gauge	Length	Tip angle
AETR18035	18G	35 mm	16°
AETR18050	18G	50 mm	16°
AETR18070	18G	70 mm	16°
AETR18090	18G	90 mm	16°
AETR18100	18G	100 mm	16°
AETR18120	18G	120 mm	16°
AETR18150	18G	150 mm	16°
AETR21090	21G	90 mm	16°
AETR21120	21G	120 mm	16°
AETR21150	21G	150 mm	16°
AETR22050	22G	50 mm	16°
AETR22070	22G	70 mm	16°
AETR22090	22G	90 mm	16°
AETR22100	22G	100 mm	16°
AETR22120	22G	120 mm	16°
AETR22150	22G	150 mm	16°
AETR23035	23G	35 mm	16°
AETR23070	23G	70 mm	16°
AETR23100	23G	100 mm	16°
AETR24035	24G	35 mm	16°

Stimy Ultra 20 - pp.07

CODE	Gauge	Length	Tip angle
AEYR18035	18G	35 mm	20°
AEYR18050	18G	50 mm	20°
AEYR18070	18G	70 mm	20°
AEYR18090	18G	90 mm	20°
AEYR18100	18G	100 mm	20°
AEYR18120	18G	120 mm	20°
AEYR18150	18G	150 mm	20°
AEYR21090	21G	90 mm	20°
AEYR21120	21G	120 mm	20°
AEYR21150	21G	150 mm	20°
AEYR22050	22G	50 mm	20°
AEYR22070	22G	70 mm	20°
AEYR22090	22G	90 mm	20°
AEYR22100	22G	100 mm	20°
AEYR22120	22G	120 mm	20°
AEYR22150	22G	150 mm	20°
AEYR23035	23G	35 mm	20°
AEYR23070	23G	70 mm	20°
AEYR23100	23G	100 mm	20°
AEYR24035	24G	35 mm	20°

Stimy Ultra 30 - pp.07

CODE	Gauge	Length	Tip angle
AEQR18035	18G	35 mm	30°
AEQR18050	18G	50 mm	30°
AEQR18070	18G	70 mm	30°
AEQR18090	18G	90 mm	30°
AEQR18100	18G	100 mm	30°
AEQR18120	18G	120 mm	30°
AEQR18150	18G	150 mm	30°
AEQR21090	21G	90 mm	30°
AEQR21120	21G	120 mm	30°
AEQR21150	21G	150 mm	30°

AEQR22050	22G	50 mm	30°
AEQR22070	22G	70 mm	30°
AEQR22090	22G	90 mm	30°
AEQR22100	22G	100 mm	30°
AEQR22120	22G	120 mm	30°
AEQR22150	22G	150 mm	30°
AEQR23035	23G	35 mm	30°
AEQR23070	23G	70 mm	30°
AEQR23100	23G	100 mm	30°
AEQR24035	24G	35 mm	30°

Stimy Plus 16° - pp.06

CODE	Gauge	Length	Tip angle
AETM18035	18G	35 mm	16°
AETM18050	18G	50 mm	16°
AETM18070	18G	70 mm	16°
AETM18090	18G	90 mm	16°
AETM18100	18G	100 mm	16°
AETM18120	18G	120 mm	16°
AETM18150	18G	150 mm	16°
AETM21090	21G	90 mm	16°
AETM21120	21G	120 mm	16°
AETM21150	21G	150 mm	16°
AETM22050	22G	50 mm	16°
AETM22070	22G	70 mm	16°
AETM22090	22G	90 mm	16°
AETM22100	22G	100 mm	16°
AETM22120	22G	120 mm	16°
AETM22150	22G	150 mm	16°
AETM23035	23G	35 mm	16°
AETM23070	23G	70 mm	16°
AETM23100	23G	100 mm	16°
AETM24035	24G	35 mm	16°

Stimy Plus 30° - pp.06

CODE	Gauge	Length	Tip angle
AEQM18035	18G	35 mm	30°
AEQM18050	18G	50 mm	30°
AEQM18070	18G	70 mm	30°
AEQM18090	18G	90 mm	30°
AEQM18100	18G	100 mm	30°
AEQM18120	18G	120 mm	30°
AEQM18150	18G	150 mm	30°
AEQM21090	21G	90 mm	30°
AEQM21120	21G	120 mm	30°
AEQM21150	21G	150 mm	30°
AEQM22050	22G	50 mm	30°
AEQM22070	22G	70 mm	30°
AEQM22090	22G	90 mm	30°
AEQM22100	22G	100 mm	30°
AEQM22120	22G	120 mm	30°
AEQM22150	22G	150 mm	30°
AEQM23035	23G	35 mm	30°
AEQM23070	23G	70 mm	30°
AEQM23100	23G	100 mm	30°
AEQM24035	24G	35 mm	30°

Stimy Echo 16° - pp.08

CODE	Gauge	Length	Tip angle
AETMN18035	18G	35 mm	16°
AETMN18050	18G	50 mm	16°
AETMN18070	18G	70 mm	16°
AETMN18090	18G	90 mm	16°
AETMN18100	18G	100 mm	16°
AETMN18120	18G	120 mm	16°
AETMN18150	18G	150 mm	16°
AETMN21090	21G	90 mm	16°
AETMN21120	21G	120 mm	16°
AETMN21150	21G	150 mm	16°
AETMN22050	22G	50 mm	16°

AETMN22050	22G	50 mm	16°
AETMN22070	22G	70 mm	16°
AETMN22090	22G	90 mm	16°
AETMN22100	22G	100 mm	16°
AETMN22120	22G	120 mm	16°
AETMN22150	22G	150 mm	16°
AETMN23035	23G	35 mm	16°
AETMN23070	23G	70 mm	16°
AETMN23100	23G	100 mm	16°
AETMN24035	24G	35 mm	16°

Stimy Echo 30° - pp.08

CODE	Gauge	Length	Tip angle
AEQMN18035	18G	35 mm	30°
AEQMN18050	18G	50 mm	30°
AEQMN18070	18G	70 mm	30°
AEQMN18090	18G	90 mm	30°
AEQMN18100	18G	100 mm	30°
AEQMN18120	18G	120 mm	30°
AEQMN18150	18G	150 mm	30°
AEQMN21090	21G	90 mm	30°
AEQMN21120	21G	120 mm	30°
AEQMN21150	21G	150 mm	30°
AEQMN22050	22G	50 mm	30°
AEQMN22070	22G	70 mm	30°
AEQMN22090	22G	90 mm	30°
AEQMN22100	22G	100 mm	30°
AEQMN22120	22G	120 mm	30°
AEQMN22150	22G	150 mm	30°
AEQMN23035	23G	35 mm	30°
AEQMN23070	23G	70 mm	30°
AEQMN23100	23G	100 mm	30°
AEQMN24035	24G	35 mm	30°

P-RF SHARP ACTIVE LENGHT 3mm - pp.10

CODE	Gauge	Length	Tip
RF1860R3S	18G	60 mm	STRAIGHT
RF1860C3S	18G	60 mm	CURVED
RF18100R3S	18G	100 mm	STRAIGHT
RF18100C3S	18G	100 mm	CURVED
RF18150R3S	18G	150 mm	STRAIGHT
RF18150C3S	18G	150 mm	CURVED
RF2060R3S	20G	60 mm	STRAIGHT
RF2060C3S	20G	60 mm	CURVED
RF20100R3S	20G	100 mm	STRAIGHT
RF20100C3S	20G	100 mm	CURVED
RF20150R3S	20G	150 mm	STRAIGHT
RF20150C3S	20G	150 mm	CURVED
RF2160R3S	21G	60 mm	STRAIGHT
RF2160C3S	21G	60 mm	CURVED
RF21100R3S	21G	100 mm	STRAIGHT
RF21100C3S	21G	100 mm	CURVED
RF21150R3S	21G	150 mm	STRAIGHT
RF21150C3S	21G	150 mm	CURVED
RF2260R3S	22G	60 mm	STRAIGHT
RF2260C3S	22G	60 mm	CURVED

P-RF SHARP ACTIVE LENGHT 5mm - pp.10

CODE	Gauge	Length	Tip
RF1860R5S	18G	60 mm	STRAIGHT
RF1860C5S	18G	60 mm	CURVED
RF18100R5S	18G	100 mm	STRAIGHT
RF18100C5S	18G	100 mm	CURVED
RF18150R5S	18G	150 mm	STRAIGHT
RF18150C5S	18G	150 mm	CURVED
RF2060R5S	20G	60 mm	STRAIGHT
RF2060C5S	20G	60 mm	CURVED
RF20100R5S	20G	100 mm	STRAIGHT
RF20100C5S	20G	100 mm	CURVED
RF20150R5S	20G	150 mm	STRAIGHT

RF20150C5S	20G	150 mm	CURVED
RF2160R5S	21G	60 mm	STRAIGHT
RF2160C5S	21G	60 mm	CURVED
RF21100R5S	21G	100 mm	STRAIGHT
RF21100C5S	21G	100 mm	CURVED
RF21150R5S	21G	150 mm	STRAIGHT
RF21150C5S	21G	150 mm	CURVED
RF2260R5S	22G	60 mm	STRAIGHT
RF2260C5S	22G	60 mm	CURVED
RF22100R5S	22G	100 mm	STRAIGHT
RF22100C5S	22G	100 mm	CURVED
RF22150R5S	22G	150 mm	STRAIGHT
RF22150C5S	22G	150 mm	CURVED

P-RF SHARP ACTIVE LENGHT 10mm - pp.10

CODE	Gauge	Length	Tip
RF1860R10S	18G	60 mm	STRAIGHT
RF1860C10S	18G	60 mm	CURVED
RF18100R10S	18G	100 mm	STRAIGHT
RF18100C10S	18G	100 mm	CURVED
RF18150R10S	18G	150 mm	STRAIGHT
RF18150C10S	18G	150 mm	CURVED
RF2060R10S	20G	60 mm	STRAIGHT
RF2060C10S	20G	60 mm	CURVED
RF20100R10S	20G	100 mm	STRAIGHT
RF20100C10S	20G	100 mm	CURVED
RF20150R10S	20G	150 mm	STRAIGHT
RF20150C10S	20G	150 mm	CURVED
RF2160R10S	21G	60 mm	STRAIGHT
RF2160C10S	21G	60 mm	CURVED
RF21100R10S	21G	100 mm	STRAIGHT
RF21100C10S	21G	100 mm	CURVED
RF21150R10S	21G	150 mm	STRAIGHT
RF21150C10S	21G	150 mm	CURVED
RF2260R10S	22G	60 mm	STRAIGHT
RF2260C10S	22G	60 mm	CURVED
RF22100R10S	22G	100 mm	STRAIGHT
RF22100C10S	22G	100 mm	CURVED
RF22150R10S	22G	150 mm	STRAIGHT
RF22150C10S	22G	150 mm	CURVED

P-RF BLUNT ACTIVE LENGHT 3mm - pp.10

CODE	Gauge	Length	Tip
RF1860R3B	18G	60 mm	STRAIGHT
RF1860C3B	18G	60 mm	CURVED
RF18100R3B	18G	100 mm	STRAIGHT
RF18100C3B	18G	100 mm	CURVED
RF18150R3B	18G	150 mm	STRAIGHT
RF18150C3B	18G	150 mm	CURVED
RF2060R3B	20G	60 mm	STRAIGHT
RF2060C3B	20G	60 mm	CURVED
RF20100R3B	20G	100 mm	STRAIGHT
RF20100C3B	20G	100 mm	CURVED
RF20150R3B	20G	150 mm	STRAIGHT
RF20150C3B	20G	150 mm	CURVED
RF2160R3B	21G	60 mm	STRAIGHT
RF2160C3B	21G	60 mm	CURVED
RF21100R3B	21G	100 mm	STRAIGHT
RF21100C3B	21G	100 mm	CURVED
RF21150R3B	21G	150 mm	STRAIGHT
RF21150C3B	21G	150 mm	CURVED
RF2260R3B	22G	60 mm	STRAIGHT
RF2260C3B	22G	60 mm	CURVED
RF22100R3B	22G	100 mm	STRAIGHT
RF22100C3B	22G	100 mm	CURVED
RF22150R3B	22G	150 mm	STRAIGHT
RF22150C3B	22G	150 mm	CURVED

P-RF BLUNT ACTIVE LENGHT 5mm - pp.10

CODE	Gauge	Length	Tip
RF1860R5B	18G	60 mm	STRAIGHT
RF1860C5B	18G	60 mm	CURVED
RF18100R5B	18G	100 mm	STRAIGHT
RF18100C5B	18G	100 mm	CURVED
RF18150R5B	18G	150 mm	STRAIGHT
RF18150C5B	18G	150 mm	CURVED
RF2060R5B	20G	60 mm	STRAIGHT
RF2060C5B	20G	60 mm	CURVED
RF20100R5B	20G	100 mm	STRAIGHT
RF20100C5B	20G	100 mm	CURVED
RF20150R5B	20G	150 mm	STRAIGHT
RF20150C5B	20G	150 mm	CURVED
RF2160R5B	21G	60 mm	STRAIGHT

RF2160C5B	21G	60 mm	CURVED
RF21100R5B	21G	100 mm	STRAIGHT
RF21100C5B	21G	100 mm	CURVED
RF21150R5B	21G	150 mm	STRAIGHT
RF21150C5B	21G	150 mm	CURVED
RF2260R5B	22G	60 mm	STRAIGHT
RF2260C5B	22G	60 mm	CURVED
RF22100R5B	22G	100 mm	STRAIGHT
RF22100C5B	22G	100 mm	CURVED
RF22150R5B	22G	150 mm	STRAIGHT
RF22150C5B	22G	150 mm	CURVED

P-RF BLUNT ACTIVE LENGHT 10mm - pp.10

CODE	Gauge	Length	Tip
RF1860R10B	18G	60 mm	STRAIGHT
RF1860C10B	18G	60 mm	CURVED
RF18100R10B	18G	100 mm	STRAIGHT
RF18100C10B	18G	100 mm	CURVED
RF18150R10B	18G	150 mm	STRAIGHT
RF18150C10B	18G	150 mm	CURVED
RF2060R10B	20G	60 mm	STRAIGHT
RF2060C10B	20G	60 mm	CURVED
RF20100R10B	20G	100 mm	STRAIGHT
RF20100C10B	20G	100 mm	CURVED
RF20150R10B	20G	150 mm	STRAIGHT
RF20150C10B	20G	150 mm	CURVED
RF2160R10B	21G	60 mm	STRAIGHT
RF2160C10B	21G	60 mm	CURVED
RF21100R10B	21G	100 mm	STRAIGHT
RF21100C10B	21G	100 mm	CURVED
RF21150R10B	21G	150 mm	STRAIGHT
RF21150C10B	21G	150 mm	CURVED
RF2260R10B	22G	60 mm	STRAIGHT
RF2260C10B	22G	60 mm	CURVED
RF22100R10B	22G	100 mm	STRAIGHT
RF22100C10B	22G	100 mm	CURVED
RF22150R10B	22G	150 mm	STRAIGHT
RF22150C10B	22G	150 mm	CURVED

