

Palium Needle Coaxial

Tru-cut needle for soft tissues biopsy to be used with automatic gun type Palium® and Magnum® Bard® with coaxial needle




Coaxial intoducer

Palium Needle suitable for using with Palium reusable gun and Bard Magnum® reusable gun.



The needle is supplied with a coaxial introducer studied to be perfectly coupled in diameter and length with the needle.

The tip of the coaxial needle is positioned at max 3 mm from the beginning of the notch of the stylet when the needle is armed.



The use of coaxial give the following **advantages** :

- In case of biopsy procedure in hard tissues that are difficult to cut, the use of coaxial permit to increase the speed of the advancement of the cannula by sliding inside it avoiding the friction of the parenchymal tissues.
- Performing multiple biopsies by inserting the needle into the coaxial that remain inside the target.
- Finding easier the exact point where to perform the biopsy thanks to the ecogenic marker.



Code needle with coaxial	Needle Gauge	Needle Length L	Configuration	Coaxial needle code	Coaxial needle gauge	Length Coaxial l
MGC1210	12G	100	NON STANDARD	PT11081	11G	81mm
MGC1213	12G	130	NON STANDARD	PT11111	11G	111mm
MGC1215	12G	150	NON STANDARD	PT11131	11G	131mm
MGC1216	12G	160	NON STANDARD	PT11141	11G	141mm
MGC1220	12G	200	NON STANDARD	PT11181	11G	181mm
MGC1410	14G	100	STANDARD	PT13082	13G	82mm
MGC1413	14G	130	NON STANDARD	PT13112	13G	112mm
MGC1416	14G	160	NON STANDARD	PT13142	13G	142mm
MGC1420	14G	200	NON STANDARD	PT13182	13G	182mm
MGC1610	16G	100	NON STANDARD	PT15083	15G	83mm
MGC1613	16G	130	NON STANDARD	PT15113	15G	113mm
MGC1615	16G	150	NON STANDARD	PT15133	15G	133mm
MGC1616	16G	160	NON STANDARD	PT15143	15G	143mm
MGC1620	16G	200	NON STANDARD	PT15183	15G	183mm
MGC1810	18G	100	STANDARD	PT17083	17G	83mm
MGC1813	18G	130	NON STANDARD	PT17113	17G	113mm
MGC1815	18G	150	STANDARD	PT17133	17G	133mm
MGC1816	18G	160	NON STANDARD	PT17143	17G	143mm
MGC1820	18G	200	STANDARD	PT17183	17G	183mm
MGC1825	18G	250	NON STANDARD	PT17233	17G	233mm
MGC1830	18G	300	NON STANDARD	PT17283	17G	283mm
MGC2010	20G	100	NON STANDARD	PT19084	19G	84mm
MGC2013	20G	130	NON STANDARD	PT19114	19G	114mm
MGC2016	20G	160	NON STANDARD	PT19144	19G	144mm
MGC2020	20G	200	NON STANDARD	PT19184	19G	184mm



M.D.L. Srl
Via Tavani 1/a 23014 Delebio (So) Italy
Tel +39 0342 682130 Fax +39 0342 691316
email : info@mdlssl.com web : www.mdlssl.com



SOFT TISSUE BIOPSY